

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:



Tested by TÜV Berlin  
FDA registered - USA  
CE marked

FUNCTIONS

- **Frame:** height-adjustable, made of steel painted with non-toxic epoxy paint
- **Seat:** triangular to hold the pelvis (seat/back at a 90° angle). It is padded and wrap-around, upholstered with strong, washable plastic
- **Bumpers:** elastic bumpers to protect the wall and protective balls on the upper ends of the rods
- **Wheels:** 80 mm fully-swivelling, rolling on ball-bearings, 2 with brakes
- **Handle bar:** adjustable in both height and depth
- **Abduction block:** adjustable to contain the pelvis

Please check off the interested products:  \_\_\_\_\_

MEASURES




- BIRILLO 1** extra small
- BIRILLO 2** small
- BIRILLO 3** medium
- BIRILLO 4** large

Frame colour:  
**blue**

Upholstery colours:  
**yellow-green-red (combined)**




Attention!  
Read carefully the user and  
maintenance manual



SIZE	AID WEIGHT	TURNING WIDTH	MAX CAP LOAD	A	B	C	D	E	F	G
1	Kg. 8,5	cm 67	Kg. 30	cm 18	cm 59	cm 50	cm 19-27	cm 36-50	cm 18	cm 22
2	Kg. 10	cm 73	Kg. 40	cm 19	cm 65	cm 51	cm 22-30	cm 40-56	cm 20	cm 26
3	Kg. 11	cm 81	Kg. 50	cm 21	cm 73	cm 59	cm 22-34	cm 42-62	cm 23	cm 30
4	Kg. 12,8	cm 89	Kg. 60	cm 28	cm 78	cm 60	cm 22-40	cm 50-76	cm 30	cm 38

ACCESSORIES

**812**  
Directional lock



**824**  
Transparent TABLE



**860**  
Concave HEADREST  
adjustable in height



Notes:

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

