

Tel. +39 0742 22927 Fax +39 0742 22637 www.ormesa.com info@ormesa.com

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:



Certified by **CERAH France**  
FDA registered - **USA**  
CE marked



### FUNCTIONS

- **Lightweight aluminium frame**, painted with non-toxic epoxy powders. It folds just like an umbrella
- **Semi-rigid backrest**, reclines in three positions with zippers
- **Semi-rigid seat bottom**, slightly inclined
- **Padded upholstery** removable for washing with removable aluminium reinforcement ribs; slots for attaching accessories
- **Net basket** on the back of the backrest
- **Removable footrest** height adjustable
- **Legrest strap**
- **Puncture-proof wheels** sliding on ball bearings: front 1.5-cm swivelling, rear 20-cm with independent brakes


Please check off the interested products:  \_\_\_\_\_

### MEASURES

- CLIP-1**
- CLIP-2**
- CLIP-3**
- Upholstery colour: **Blue and grey combined**  
Frame colour: **Grey**



**Attention!**  
Read carefully the user and maintenance manual



SIZES	MAX CAP. LOAD	AID WEIGHT	A	B	C	D	E	F	G	H	I	L	M	N	O	P	Q
1	<b>Kg. 30</b>	Kg. 10,8	cm 40	cm 30	cm 51	cm 99	0°-12°	cm 25	cm 55	18°	cm 46	cm 15-23	cm 73	cm 80	cm 110	cm 30	cm 32
2	<b>Kg. 35</b>	Kg. 11	cm 40	cm 33	cm 51	cm 99	0°-12°	cm 30	cm 59	18°	cm 48	cm 21-31	cm 73	cm 90	cm 110	cm 30	cm 32
3	<b>Kg. 40</b>	Kg. 11,2	cm 40	cm 36	cm 51	cm 104	0°-12°	cm 35	cm 64	18°	cm 49	cm 27-37	cm 79	cm 96	cm 115	cm 32	cm 35

### ACCESSORIES

**812**  
Wheel direction locks

**816**  
Side supports for the trunk

**817**  
Leg straps

**818**  
Thermal cover

**819** Canopy

**825**  
Rain cover (it can only be attached after the canopy 819)

**827**  
Foot straps

**828**  
Chest belt

**835**  
Adjustable headrest

**853**  
Vest harness

**891** Set of 4 tie down hooks for the forward-facing transport of the pushchair in a motor vehicle (private vans/cars, buses..)

**894**  
45° Pelvic belt

**903**  
Five point vest harness

**906**  
Five point harness

Notes:

.....

.....

.....

.....



Congratulations for choosing **ORMESA!**