

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:

**F**  
**U**  
**N**  
**C**  
**T**  
**I**  
**O**  
**N**  
**S**

- Steel and aluminum folding **frame** with a non-toxic epoxy varnish finish
- Fire-retardant padding and removable **upholstery**
- Rigid **bakrest** with 4 quick release recline positions
- Fixed rigid **seat**
- Adjustable and tip-up **footrest**
- Drum **brakes**
- **Wheels:** front puncture proof cm 17,5 castor wheels, and back pneumatic cm 30 wheels, acting as dampers.




Certified by **CERAH France**  
€marked

Please check off the interested products:

**M**  
**E**  
**A**  
**S**  
**U**  
**R**  
**E**  
**S**

- OBI-2**
- OBI-3**
- OBI-4**

Frame colour: **blue**  
Upholstery colour: **green-pink-blue** (combined)



SIZE	max load cap	weight of the aid	A	B	C	D	E	G	H	I	L	M
2	Kg. 50	Kg. 17,4	cm 62	cm 33	cm 30	cm 23-35	18°	cm 60	cm 97	cm 99	cm 92	cm 50
3	Kg. 60	Kg. 18,5	cm 68	cm 37	cm 36	cm 29-44	18°	cm 62,5	cm 103	cm 109	cm 95	cm 46
4	Kg. 70	Kg. 20,1	cm 76	cm 41	cm 40	cm 32-47	18°	cm 66,5	cm 105	cm 124	cm 103	cm 51

**A**  
**C**  
**C**  
**E**  
**S**  
**S**  
**O**  
**R**  
**I**  
**E**  
**S**

<input type="checkbox"/> <b>812</b> Wheel direction locks 	<input type="checkbox"/> <b>816</b> Side supports 	<input type="checkbox"/> <b>817</b> Leg straps 	<input type="checkbox"/> <b>818</b> Thermal cover 	<input type="checkbox"/> <b>819</b> Canopy 	<input type="checkbox"/> <b>824</b> Table with edge and hollow 
<input type="checkbox"/> <b>825</b> Rain cover (it can only be attached after the canopy 819) 	<input type="checkbox"/> <b>827</b> Foot straps 	<input type="checkbox"/> <b>828</b> Seat belt 	<input type="checkbox"/> <b>834</b> Abduction block 	<input type="checkbox"/> <b>835</b> Headrest with two side elements 	<input type="checkbox"/> <b>837</b> Tray 
<input type="checkbox"/> <b>839</b> Front handle 	<input type="checkbox"/> <b>853</b> Adjustable harness 	<input type="checkbox"/> <b>858</b> Basket 	<input type="checkbox"/> <b>891</b> Set of 4 tie down hooks for the forward-facing transport of the pushchair in a motor vehicle (private vans/cars, buses...)  <b>CRASH TESTED</b> according to the ISO 7176/19	<input type="checkbox"/> <b>894</b> 45° Pelvic belt 	

Notes: .....

