

INGLESE - cod.220491 - 23/09/2022



ORMESA®
• MADE IN ITALY SINCE 1980 •

multifunction
Table



multifunction Table

Perfect Stability.

The **Multifunction Table** by Ormesa has been designed for **training and work activities and as a complement to other devices such as chairs and/or baby chairs.**

It provides a wide table top, ergonomic and functional to exercise manual activities, reading, allowing the user to stabilize his/her trunk in order to facilitate the functional use of the upper limbs, arrange a stable posture that allows a wide and structured view of the space, distances and environment surrounding the subject in his sitting position.

The **Multifunction Table is available in three different sizes:** the small size for children, the medium size for teenagers and the large size for adults.

PRODUCT
MARKED 

MULTIFUNCTION TABLE is a MEDICAL DEVICE that COMPLIES WITH THE SAFETY REQUIREMENTS defined with the REGULATION (EU) 2017/745



Maximum convenience

Transport and storage at your fingertips

New design

Light and foldable for transport and storage. Its low weight allows you **to carry even the largest size as if it were a briefcase.**

It has a wide base to allow easy access of the wheelchair/baby chair. The feet are built with non-slip material for more safety.

Frame painted with non-toxic epoxy powders, after antioxidant treatment.



Functional top

Ormesa design Table

Ergonomic

Table top with large surface, suitably shaped with recess to support the forearms and facilitate the approach of the person.

Made of a **highly resistant material, easy to sanitize and hygienic**. Design without edges that does not trap food residues.

Adjustable in tilt with safety **system that prevents the frontal overturning of the top**.



Easily
sanitizable





Technical features



Height adjustment system

Simple, reliable and precise adjustment system: the levers allow you to easily tighten the desired height while the internal friction system prevents the accidental fall of the table top.

Tools free

Adjustments without the use of tools but thanks to levers, safe and easy to operate.



Bookend and pen holder

Also useful with the inclined table top.




Hook

Lateral Hook to store backpack or bag.



Table sizes and weights


	small	medium	large
A: table width	64 cm ~ 25,1 in	80 cm ~ 31,4 in	80 cm ~ 31,4 in
B: base width	63 cm ~ 24,8 in	80 cm ~ 31,4 in	80 cm ~ 31,4 in
C: overall length	54 cm ~ 21,2 in	61 cm ~ 24 in	61 cm ~ 24 in
D: ground height - table top	from 44 to 57 cm ~ from 17,3 to 22,4 in	from 53 to 77 cm ~ from 20,8 to 30,3 in	from 69 to 108 cm ~ from 27,1 to 42,5 in
E: table height when closed	17 cm ~ 6,6 in	16 cm ~ 6,2 in	16 cm ~ 6,2 in
F: table length when closed	84 cm ~ 33 in	104 cm ~ 40,9 in	120 cm ~ 47,2 in
G: table tilt (recommended)	from 0° to 60°	from 0° to 60°	from 0° to 60°
I: recess width	38,3 cm ~ 15 in	42,4 cm ~ 16,6 in	42,4 cm ~ 16,6 in
L: recess depth	16 cm ~ 6,2 in	14 cm ~ 5,5 in	14 cm ~ 5,5 in
H: user height	from 75 to 90 cm ~ from 29,5 to 35,4 in	from 85 to 120 cm ~ from 33,4 to 47,2 in	from 115 to 170 cm ~ from 45,2 to 66,9 in
Total Weight	6,7 kg ~ 14,7 lb	8,2 kg ~ 18 lb	9,2 kg ~ 20,2 lb











CE complies with MDR (EU) 2017/745

COMPANY WITH
QUALITY SYSTEM
CERTIFIED BY DNV GL
= **ISO 13485** =

Ormesa srl - Via delle Industrie 6/8 | 06034 Foligno (PG) - ITALY
Tel. +39 0742 22927 | FAX +39 0742 22637 | info@ormesa.com

www.ormesa.com