



ORMESA®
• MADE IN ITALY SINCE 1980 •

NOVUS



NOVUS

Maximum seating comfort.

Novus by Ormesa is a **stroller for children and teens** with postural needs **available in three sizes.**

Equipped with a **folding aluminum frame**, it has been designed to ensure **maximum comfort** to the user: the seat is slightly tilted, while the cushion, padded and removable, is **highly breathable**. The **wheels are provided with drum brakes**; they can be easily removed with a snap-out system, to facilitate the transport.

Thanks to a wide range of components, Novus adapts to any needs.



Novus passed **CRASH TEST** according to **ANSI RESNA WC/19**

PRODUCT MARKED **CE**

NOVUS is a **MEDICAL DEVICE** that **COMPLIES WITH THE SAFETY REQUIREMENTS** defined with the **REGULATION (EU) 2017/745**




Technical features and adjustments



Foldable frame

made in aluminum. Painted with non-toxic poxy powders.



Removable upholstery

padded, ergonomic and fixed to the seat thanks to press-studs. It was designed to make the change of position of the side supports easier, and to fix or unfix the accessories without causing any trouble to the child while seated. It grows up with the frame.



CERTIFICATE
MED annex A.1/3.19-20-21



Backrest

adjustable in height and tilt.

Seat

adjustable in depth and width thanks to the padded sides to the pelvis.



Legrest

adjustable in tilt.

Flip-up footrest

adjustable in height.



Removable wheels

removable, sliding on ball bearings.

Front wheels: 17,5 cm, swivel and puncture-proof.

Rear wheels: 30 cm, pneumatic with **independent drum brakes**.
Puncture-proof rubber on request.



Components



852 Headrest with parietal supports

adjustable in height and width.



863 Headrest with occipital parietal supports

adjustable in height, width, tilt and forward-backward.



853 Vest harness



903 Five point vest harness



906 Five point harness



838 Lateral supports
adjustable in height and width.

Components



868 Wrappable and flexible trunk supports
Adjustable in height and width.



894 45° pelvic belt



844 Side protecting covers



834 Abduction block
adjustable in depth.



834-R Abduction block
Adjustable in depth and abduction.



824 Transparent tray
with recess.

Components



839 Front handle



827 Foot straps



892 Covering for footrest



812 Directional locks for front wheels



818 Thermal cover




819 Canopy

Components



825 Rain cover
to be fixed to the 819 canopy.



905 Hand brake




922 Bottle holder



858 Net basket for 869 base



911 Vent tray




891 Set of 4 tie down hooks
for the forward-facing transport of
the pushchair in a motor vehicle
(private cars, buses..).


*See User and Maintenance Manual
provided with the aid.*

Novus sizes and weight

novus 4

A: backrest width	46 cm	18,1 in
B: seat bottom width	from 28 to 38 cm	from 11 to 14,9 in
C: overall width	66 cm	25,9 in
D: overall height	107 cm	42,1 in
E: overall length	96 cm	37,7 in
F: seat bottom depth	from 39 to 45 cm	from 15,3 to 17,7 in
G: backrest height	from 75 to 81 cm	from 29,5 to 31,8 in
I: seat to footrest height	from 30 to 44 cm	from 11,8 to 17,3 in
L: folded stroller length without wheels	110 cm	43,3 in
M: folded stroller height without wheels	42 cm	16,5 in
N: folded stroller length with wheels	90 cm	35,4 in
O: folded stroller height with wheels	56 cm	22 in
H: user height	from 145 to 170 cm	from 57 to 66,9 in
Total weight with wheels	22,1 kg	48,7 lb
Total weight without wheels	18,5 kg	40,7 lb
Max capacity	65 kg	143,3 lb





 complies with MDR (EU) 2017/745

COMPANY WITH
QUALITY SYSTEM
CERTIFIED BY DNV GL
= ISO 13485 =

Ormesa srl - Via delle Industrie 6/8 | 06034 Foligno (PG) - ITALY
Tel. +39 0742 22927 | FAX +39 0742 22637 | info@ormesa.com

www.ormesa.com