

DISTRIBUTOR/AGENT	DATE OF ORDER
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:



ISO 7176/19 and  
ANSI RESNA WC/19



very breathable  
and fire retardant  
upholstery



Marked



Attention!  
Read carefully the user and  
maintenance manual

### FEATURES

- **More sizes in one frame**
- **Frame:** light and umbrella folding system, in aluminium. It is painted with non-toxic epoxy powders. When folded Trollino is compact and it allows you to carry it as if it was a travel trolley.
- **Push handles with manual and parking brakes.**
- Resistant and durable **backrest** guaranteed by the cross on the back. Reclinable on three positions and easily adjustable in height.
- Slightly tilted **seat**. It is easily adjustable in depth and width thanks to padded pelvic side supports (additional components 958). The cross below ensures stability, comfort and durability.
- Extremely breathable **upholstery** (not removable). It is resistant and easily washable. **Fire Retardant** according to *EN 1021-1:2014 and EN 1021-2:2014*.
- User manual **bag** and wrench on the back of the backrest.
- **Legrest** adjustable in height.
- Elevatable and swivel **footrests**, without encumbrance.
- **Solid wheels:** 17 cm (7 in) front swiveling wheels and 25 cm (10 in) rear wheels with **independent drum brakes**. **Overcoming step bar** on both rear wheels.

### Please check off the interested products




Backrest colour:      Seat colour: **Black**      Frame colour: **Black**

- Red
- Pink
- Blue
- Black

### MEASURES

Unit of Measurement: cm (in)



**A: 36 cm (14 in).** For different width (32 and 28 cm - 12,6 and 11 in) please select 958 padded pelvic side supports.  
\*Height adjustment range with add-on 951 max. 105 cm (41,3 in).

C  
O  
M  
P  
O  
N  
E  
N  
T  
S

**951** Height adjustable push handles. Tools free!  
**COMPONENT TO BE ORDERED WITH THE PRODUCT**

↙ **COMPONENTS TO CHOOSE ALSO AFTERWARDS**

**835** Shaped and enveloping headrest adjustable in height.  
 S-M

**865** Padded headrest adjustable in height.  
 S-M

**838** Trunk side supports adjustable in height, width and inclination.  
 mini  S

**868** Wrappable and flexible trunk support, multiadjustable in height, width and inclination.  
 S

**958** Padded pelvic side supports 32 cm (12,6 in) seat width.  
 **Extra Padding** 28 cm (11 in) seat width

**947** Pelvic belt with variable angle.  
 S  M

**920** 4 point pelvic belt  
 S  M

**853** Vest harness.  
 mini  S

**853 slim** 4 point shaped harness.  
 M

**834N** Padded and removable abduction block.

**839** Front handlebar adjustable in tilt.

**839** Front handlebar adjustable in tilt.

**960** Heel rests.

**956** Footrests adjustable in tilt.

**858** Shopping basket.

**819** Elevating Canopy with bag. **Foldable with the frame.** **NEW**

**825** Rain cover to fix on the 819 canopy. **NEW**

**825B** Back and lateral protection clothe to fix on the 819 canopy. **NEW**

**818** Thermic cover.

**968** Padded covers for armrests. Easily removable for washing. **NEW**

**844** Padded side protection covers.

**812** Front wheel direction locks.

**891** Set of 4 tie down hooks for the forward-facing transport of the pushchair in a motor vehicle. **CRASH TESTED**

COMPONENTS MEASURES	<b>947</b>	<b>920</b>	<b>838</b>	
	<p><b>small</b> A=19 cm/7,4 in B=6 cm/2,3 in</p> <p><b>medium</b> A=25,5 cm/10 in B=8 cm/3,1 in</p>	<p><b>small</b> A=16 cm/6,2 in B=15,5 cm/6,1 in</p> <p><b>medium</b> A=23 cm/9 in B=19 cm/7,4 in</p>	<p><b>mini</b> A=9 cm/3,5 in B=7 cm/2,7 in</p> <p><b>small</b> A=10 cm/3,9 in B=12 cm/4,7 in</p>	C= 14-23 cm/5,5-9 in
<b>868</b>	<b>853</b>	<b>853 slim</b>	<b>834N</b>	
<p><b>small</b> A=26 cm/10,2 in B=9 cm/3,5 in</p>	<p><b>mini</b> A=20,5 cm/8 in B=28 cm/11 in</p> <p><b>small</b> A=25 cm/9,8 in B=34 cm/13,3 in</p>	<p><b>medium</b> A=33 cm/12,9 in B=46 cm/18,1 in C=5 cm/1,9 in</p>	<p><b>mini/S</b> A=12 cm/4,7 in B=4 cm/1,5 in C=14 cm/5,5 in</p>	
<b>951</b>				

Notes.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....