

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:



according to the ISO 7176/19



### FUNCTIONIONS

- Lightweight aluminium **frame**, painted with non-toxic epoxy powders. It **folds just like an umbrella**.
- **Semi-rigid backrest**, reclines in three positions with zippers.
- **Height adjustable Handles, Tools free**.
- **Semi-rigid seat** bottom, slightly inclined.
- **Padded upholstery** removable for washing with **removable aluminium reinforcement ribs**. Extremely breathable and quick to dry, excellent for **load distribution**; **Fireproof fabric** according to EN 1021-1:2014 and EN 1021-2:2014.
- **Net basket** on the back of the backrest.
- Removable **footrest** height adjustable.
- **Legrest strap**.
- Puncture-proof **wheels**: front 15 cm swivelling, rear 20 cm with independent brakes **and overcoming step bar**.

Attention! Read carefully the user and maintenance manual


### Please check off the interested products:



Upholstery colour: **Black and red** Frame colour: **Black**

- CLIP 1**
- CLIP 2**
- CLIP 3**

Measures in cm (Inches):



SIZE	Max. Load	Total Weight	A	B	C	D	E	F	G	H	I	L	M	N	O	P	Q
1	30 kg 66,1 lb	10,3 kg 22,7 lb	40 (15,7)	51 (20)	101-108 (39,7-42,5)	0°-12°	25 (9,8)	55 (21,6)	18°	46 (18,1)	12-20 (4,7-7,8)	78 (30,7)	84 (33)	114 (44,8)	30 (11,8)	34 (13,3)	110-125 (43,3-49,2)
2	35 kg 77,1 lb	10,5 kg 23,1 lb	40 (15,7)	51 (20)	101-108 (39,7-42,5)	0°-12°	30 (11,8)	59 (23,2)	18°	48 (18,1)	21-31 (8,2-12,2)	78 (30,7)	92 (36,2)	114 (44,8)	30 (11,8)	34 (13,3)	120-140 (47,2-55,1)
3	40 kg 88,1 lb	10,7 kg 23,5 lb	40 (15,7)	51 (20)	106-113 (41,7-44,4)	0°-12°	35 (13,7)	63 24,8	18°	49 (19,2)	25-35 (9,8-13,7)	84 (33)	98 (38,5)	118 (46,4)	34 (13,3)	36 (14)	135-155 (53,1-61)

ADDITIONAL COMPONENTS

**835** Headrest with parietal supports.



**816** Trunk lateral containment pads.



**853** Vest harness. (S for CLIP S, M, M for CLIP L)



**906** Five point harness. (S for CLIP S, M, M for CLIP L)



**903** Five point vest harness. (S for CLIP S, M, M for CLIP L)



**894** 45° Pelvic belt. (S for CLIP S, M, M for CLIP L)



**834N** Padded and removable abduction block.  mini/S  M



**858** Shopping basket.



**827** Foot straps.



**812** Front wheel direction locks.



**819** Canopy.



**818** Thermal cover.




**825** Rain cover. (to fix to the canopy 819)



**891** Set of 4 tie down hooks for the forward-facing transport of the pushchair in a motor vehicle (private vans/cars, buses..)




**853**



Small A=30 cm/11,8 in B=31,5 cm/12,4 in	Medium A=39,5 cm/15,5 in B=39 cm/15,3 in
---	--

**834N**



mini/S A=12 cm/4,7 in B=4 cm/1,5 in C=14 cm/5,5 in	Medium A=15 cm/5,9 in B=4 cm/1,5 in C=17 cm/6,6 in
---	---

Notes.....

.....

.....

.....

.....