

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:

FUNCTIONS

- Lightweight aluminium **frame**, painted with non-toxic epoxy powders. It folds just like an umbrella.
- **Semi-rigid backrest**, reclines in three positions with zippers.
- **Height adjustable Handles, Tools free.**
- **Semi-rigid seat** bottom, slightly inclined.
- **Padded upholstery** removable for washing with **removable aluminium reinforcement ribs**. Extremely breathable and quick to dry, excellent for **load distribution**; **Fireproof fabric** according to EN 1021-1:2014 and EN 1021-2:2014.
- **Net basket** on the back of the backrest.
- Removable **footrest** height adjustable and foldable.
- **Legrest strap.**
- Puncture-proof **wheels**: front 17,5 cm swivelling, rear 25 cm with independent brakes **and overcoming step bar.**



CE marked

Attention! Read carefully the user and maintenance manual

Please check off the interested products:

Upholstery colour: **Black and grey** Frame colour: **Black**

CLIP 4

Measures in cm (Inches):

MEASURES



Max. load	Total weight	A	B*	C	D	E	F	G	H	I	L	M	N	O	P	Q	R
76 kg 167,5 lb	13,9 kg 30,6 lb	45 (17,7)	33 (12,9)	57 (22,4)	113 (44,4)	0°-10°	38 (14,9)	71 (27,9)	20°	59 (23,2)	37-47 (14,5-18,5)	89 (35)	108 (42,5)	117 (46)	60 (23,6)	37 (14,5)	150-170 (59-66,9)

* Seat width without padded sides (B) cm 45 (17,7 in)

ADDITIONAL COMPONENTS

835 Headrest with parietal supports.



816 Trunk lateral containment pads.



853 Vest harness. (size M)



906 Five point harness. (size M)



903 Five point vest harness. (size M)



894 45° Pelvic belt. (size M)



834N Padded and removable abduction block. (size M)



858 Shopping basket.



827 Foot straps.



812 Front wheel direction locks.



819 Canopy.




818 Thermal cover.



825 Rain cover. (to fix to the canopy 819).




853



Medium
A=39,5 cm/15,5 in
B=39 cm/15,3 in

834N



Medium
A=15 cm/5,9 in
B=4 cm/1,5 in
C=17 cm/6,6 in

Notes.....
.....
.....
.....