

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:



**F
E
A
T
U
R
E
S**

- Aluminium **frame** with a non-toxic epoxy varnish finish.
- Washable special fireproof, visco-elastic plastic **covering**. It allows an excellent perspiration and it's highly wear-resistant. The sufficiently elastic covering adapts to the shape of the shoulders and of the bottom.
- Padded **headrest** adjustable in tilt, height and in forward-backward.
- Retentive **backrest**, it can be fixed in 2 positions in order to have the seat adjustable in depth (2 different ones). Reclining back thanks to the back rack locking.
- **Seat** with a 15° tilt, in order to prevent the user from slipping forward.
- Detachable and turning **armrests**, in order to facilitate the user's transfers.
- **Legrest** adjustable in tilt.
- One only overturnable **footrest**, in order to help the approach of the care giver, the sitting and the getting down of the user.
- Unpiercerable **wheels** 12 cm, back wheels with brake.



Attention!
 Read carefully the user and maintenance manual



Please check off the interested products: _____

**M
E
A
S
U
R
E
S**

203-5 SHOWER CHAIR

Frame colour: **White**
 PVC-string colour: **Crystal**

TOTAL WEIGHT KG 15
MAX LOAD KG 110



ADDITIONAL COMPONENTS

894 45° Pelvic belt

834 Abduction block adjustable in depth



Notes:

INGLESE - Cod 22589 - 30.07.21



download and send the form to info@ormesa.com