



Product registered **FDA (USA)**

Product marked **CE**

- **ALUMINUM FRAME**, foldable for easy transportation and convertible, from **front configuration to rear configuration and vice versa**.
- **PADDED PELVIC and THORACIC SUPPORTS**, independently adjustable in height, width and circumference; allow an easy access thanks to the back opening. They hold the user just like two hands and this gives him/her a sense of security and stability allowing the user to concentrate only on walking without any fear.
- **GRADUATED BARS** as a reference on height adjustments.
- **ERGONOMIC HARNESS**, ergonomic, comfortable, adjustable and easily removable for washing. It supports the users in case of necessity without generating painful pressures.
- **WHEELS**, light and puncture-proof: 2 front swiveling, Ø15 and 2 rear solid, Ø cm 20, for indoor and outdoor use.
- **DIRECTIONAL LOCKS** on front wheels.
- **ANTI-ROLLBACK** device, **ADJUSTABLE BRAKES**, parking brake on rear wheels, easy to operate by the red pedal.
- **SEMI-ELASTIC BUMPERS** on front wheels to protect against bumps.




Attention!
Read carefully the user and
maintenance manual

FRONT DRIVE

Inner circumference pelvic and thoracic support

from 53 to 79 cm
from 20,8 to 31,1 in




WEIGHT kg 15/33 lb

MAX. LOAD kg 35/77,1 lb

REAR DRIVE

Inner circumference pelvic and thoracic support

from 53 to 79 cm
from 20,8 to 31,1 in




WEIGHT kg 15/33 lb

MAX. LOAD kg 35/77,1 lb

USER DIMENSIONS:

USER HEIGHT (C):
min 75 cm/29,5 in
max 95 cm/37,4 in



IMPORTANT:

Make sure that measures A and B indicated in the figure do not exceed the suggested range indicated for the supports.



download and send the form to
info@ormesa.com



Front Drive

Rear Drive

CHOOSE THE CONFIGURATION >	Front Drive				Rear Drive		
	PT	PA	P	PTM	PT	PA	P
ADDITIONAL COMPONENTS v Including: > thoracic support > pelvic support > ergonomic harness > handlebar	<input type="checkbox"/> Including: > thoracic support > pelvic support > ergonomic harness > handlebar	<input type="checkbox"/> Including: > pelvic support > ergonomic harness > arm supports	<input type="checkbox"/> Including: > pelvic support > ergonomic harness > handlebar	<input type="checkbox"/> Including: > transparent tray > handlebar > pelvic support > ergonomic harness	<input type="checkbox"/> Including: > thoracic support > pelvic support > ergonomic harness > handles	<input type="checkbox"/> Including: > pelvic support > ergonomic harness > arm supports	<input type="checkbox"/> Including: > pelvic support > ergonomic harness > handles
865 Multiaadjustable headrest	NA	NA	NA	NA	<input type="checkbox"/>	<input type="checkbox"/> available in combination with thoracic support 925	<input type="checkbox"/> available in combination with thoracic support 925
925 Thoracic support	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> when transparent tray is removed	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
924 Pelvic support	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
930 Ergonomic harness	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
894 Rear table belt for grillo mini PTM	NA	NA	NA	<input type="checkbox"/>	NA	NA	NA
926 Adjustable handlebar	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	NA	NA	NA	NA
809 mini Arm supports	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	NA	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
809 mini Arm supports with curved bars	NA	<input checked="" type="checkbox"/>	NA	<input type="checkbox"/> when handlebar and tray are removed	NA	<input checked="" type="checkbox"/>	NA
943 mini Arm straps for size mini	<input type="checkbox"/> available in combination with arm supports 809	<input type="checkbox"/>	<input type="checkbox"/> available in combination with arm supports 809	<input type="checkbox"/> available in combination with arm supports 809	<input type="checkbox"/> available in combination with arm supports 809	<input type="checkbox"/>	<input type="checkbox"/> available in combination with arm supports 809
937 Ergonomic handles for size mini	NA	NA	NA	<input type="checkbox"/> when handlebar is removed	NA	NA	NA
927 Removable knobs	NA	NA	NA	NA	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
946 Ergonomic saddle	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
890-SC Distal abductor	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
923 Additional assistant push handle	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>