

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:



MARKED 

FEATURES

- **No edge frame**, in order to enable trunk rotation and side-arm motion. It's adjustable in height by knobs, without using wrenches, "tools free".
- **Birch wood table**. It is removable and adjustable in height, tilt and depth.
- **Trunk Support** adjustable in its circumference thanks to Velcro straps, in height, back and forward.
- **Pelvic Support** adjustable in its circumference thanks to Velcro straps, in height, back and forward.
- **Numbered bars** to find the right setting, determined by the therapist.
- **Knee-pads** adjustable in height, width, depth and in a rotatory movement.
- **Heelrests** adjustable in depth and width.
- **Birch wood footboard**. It gives a safe guide to the child.
- **Twinned solid Wheels**, 75 mm. diam., rear wheels with brakes. Thanks to the wheels the aid can be moved from one room to another having the child on it.



Attention!
 Read carefully the user and maintenance manual


Please check off the interested products:

MINI STANDY-1

MINI STANDY-2

MINI STANDY-3

MEASURES




size	Max. Load	Tot. Weight	A	B	C	D*	E	F	G	L	M**	N
1	35 kg 77,1 lb	19 kg 41,8 lb	60 cm 23,6 in	74-86 cm 29,1-33,8 in	49-77 cm 19,2-30,3 in	64 cm 25,1 in	14-26 cm 5,5-10,2 in	41-55 cm 16,1-21,6 in	53-72 cm 20,8-28,3 in	60 cm 23,6 in	18 cm 7 in	75-100 cm 29,5-39,3 in
2	45 kg 99,2 lb	21 kg 46,2 lb	70 cm 27,5 in	78-90 cm 30,7-35,4 in	59-87 cm 23,2-34,2 in	65 cm 25,5 in	14-32 cm 5,5-12,5 in	52-68 cm 20,4-26,7 in	67-94 cm 26,3-37 in	65 cm 25,5 in	22 cm 8,6 in	95-125 cm 37,4-49,2 in
3	55 kg 121,2 lb	23 kg 50,7 lb	80 cm 31,4 in	83-95 cm 32,6-37,4 in	69-100 cm 27,1-39,3 in	70 cm 27,5 in	14-42 cm 5,5-16,5 in	62-86 cm 24,4-33,8 in	82-110 cm 32,2-43,3 in	70 cm 27,5 in	29 cm 11,4 in	120-150 cm 47,2-59 in

ADDITIONAL COMPONENTS

944 Circumference reduction parts
 Thickness: 3 cm. each side



865 Multiadjustable headrest



download and send the form to info@ormesa.com

Notes: _____
