

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:



Marked

**F  
U  
N  
C  
T  
I  
O  
N  
S**

- Light and foldable **frame** for transport and storage.
- Ergonomic **table top**, with large surface and recess. Made of a highly resistant material, easy to sanitize and hygienic.
- Simple, reliable and precise **height adjustment system**, with internal **friction system** that prevents the accidental fall of the table top.
- **Tools free**, adjustments without the use of tools.
- Provided of **bookend, pen holder and lateral hook**.
- It is available in **three sizes** for children, teenagers and adults.




**Attention!**  
Read carefully the user and maintenance manual

Please check off the interested products:  \_\_\_\_\_

**M  
E  
A  
S  
U  
R  
E  
S**

- Multifunction Table **size Small**
- Multifunction Table **size Medium**
- Multifunction Table **size Large**

Frame colour: **Pebble grey**  
Table top colour: **Silk grey**



	A	B	C	D	E	F	G	I	L	User height	Total weight
<b>Small</b>	64 cm 25,1 in	63 cm 24,8 in	54 cm 21,2 in	from 44 to 57 cm from 17,3 to 22,4 in	17 cm 6,6 in	84 cm 33 in	from 0° to 60°	38,3 cm 15 in	16 cm 6,2 in	from 75 to 90 cm from 29,5 to 35,4 in	6,7 kg 14,7 lb
<b>Medium</b>	80 cm 31,4 in	80 cm 31,4 in	61 cm 24 in	from 53 to 77 cm from 20,8 to 30,3 in	16 cm 6,2 in	104 cm 40,9 in	from 0° to 60°	42,4 cm 16,6 in	14 cm 5,5 in	from 85 to 120 cm from 33,4 to 47,2 in	8,2 kg 18 lb
<b>Large</b>	80 cm 31,4 in	80 cm 31,4 in	61 cm 24 in	from 69 to 108 cm from 27,1 to 42,5 in	16 cm 6,2 in	120 cm 47,2 in	from 0° to 60°	42,4 cm 16,6 in	14 cm 5,5 in	from 115 to 170 cm from 45,2 to 66,9 in	9,2 kg 20,2 lb

Notes:  
-----  
-----  
-----

ENGLISH - cod.220495 - 05.2022

download and send the form to [info@ormesa.com](mailto:info@ormesa.com)