

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:



Certificate CERAH France  
 Marked 

### FEATURES

- Steel folding **frame**. Painted with non-toxic epoxy vanish.
- Removable **upholstery** covered with fire retardant fabric.
- **Rigid backrest** adjustable in tilt on 4 positions, thanks to a simple quick release mechanism.
- Rigid **seat**.
- **Leg rest** adjustable in tilt with quick lock mechanism.
- Tip-up **footrest** with heel pad.
- **Drum brakes**.
- **Wheels**: front wheels 17,5 cm swivel and unpierceable and rear wheels 30 cm unpierceable.

### Please check off the interested products:





MEASURES

- OBI-2  
 OBI-3  
 OBI-4

Frame Colour:  
**blue**

Upholstery Colour:  
**green-pink-blue  
 (combined)**



 **Attention!**  
 Read carefully the user and maintenance manual

SIZE	Max. Load	Total Weight	A	B	C	D	E	F	G	L	M	H
2	Kg. 50	Kg. 17,4	cm 60	cm 33	cm 99	cm 97	cm 30	cm 62	cm 23-35	cm 92	cm 50	cm 120-140
3	Kg. 60	Kg. 18,5	cm 62,5	cm 37	cm 109	cm 103	cm 36	cm 68	cm 29-44	cm 95	cm 46	cm 135-155
4	Kg. 70	Kg. 20,1	cm 66,5	cm 41	cm 124	cm 105	cm 40	cm 76	cm 32-47	cm 103	cm 51	cm 150-175

ADDITIONAL COMPONENTS

<input type="checkbox"/> <b>835</b> Headrest with parietal support	<input type="checkbox"/> <b>816</b> Side supports	<input type="checkbox"/> <b>834</b> Abduction block	<input type="checkbox"/> <b>853</b> Adjustable harness	<input type="checkbox"/> <b>894</b> 45° pelvic belt
<input type="checkbox"/> <b>827</b> Foot straps	<input type="checkbox"/> <b>839</b> Front handle	<input type="checkbox"/> <b>837</b> Tray without hollow	<input type="checkbox"/> <b>819</b> Canopy	<input type="checkbox"/> <b>825</b> Rain cover (to fix at the canopy 819)
<input type="checkbox"/> <b>818</b> Thermal cover	<input type="checkbox"/> <b>858</b> Net basket	<input type="checkbox"/> <b>812</b> Directional locks for front wheels	<input type="checkbox"/> <b>891</b> Set of 4 tie down hooks for the forward-facing transport of the pushchair in a motor vehicle (private cars, buses..).	notes: ----- ----- -----

INGLESE - Cod 22580 - 08.2021