

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:



### FEATURES

- **Adjustable handlebar** in height, backward-forward and reversible.
- **Wide and ergonomic handles**, to increase comfort.
- **Steering limiters with springs**, that allow the handlebar to return to the central position. Releasing to allow more steering. It's possible to choose between two different elastic strengths.
- **Memory foam padded trunk support**, adjustable in height, backward-forward and free-angle tilt, with breathable padded upholstery, easily removable for washing.
- **Front brake with locking device** for ascent and descent, that can be operated by the lever either on the right or on the left side of the handlebar. The distance between the brake lever and the handlebar is adjustable by screw.
- **Ergonomic saddle** adjustable in height, tilt and anterior-posterior direction.
- **Pneumatic tires** with reflector stripes.  
**Quick release front wheel.**  
**Rear wheels** easily removable with key.

Product marked



**Attention!**  
Read carefully the user and maintenance manual


Please check off the interested products:  \_\_\_\_\_

size M

size L

measures in cm and (inches)

### MEASURES






	A	B	C	D <sup>1</sup>	E	F <sup>2</sup>	G	H	Weight <sup>3</sup>	Max Load
M	55 (21,6)	86 (33,8)	84-97 (33-38,1)	159 (62,6)	53-64 (20,8-25,1)	58-73 (22,8-28,7)	106 (41,7)	120-150 (47,2-59)	17 kg (37,4 lbs)	70 kg (154 lbs)
L	62 (24,4)	92 (36,2)	95-107 (37,4-42,1)	187 (73,6)	71-84 (27,9-33)	70-92 (27,5-36,2)	127 (50)	150-190 (59-74,8)	18 kg (39,7 lbs)	100 kg (220 lbs)

<sup>1</sup> Size M with 26" wheel, D<sub>max</sub> = 164 cm (64,5") - Size L with 28" wheel, D<sub>max</sub> = 192 cm (75,5").

<sup>2</sup> With standard saddle. It is possible to further adjust the height position of the saddle (±3 cm - 1,18") by acting on the fixing clamp.

<sup>3</sup> Weight calculated without saddle abatement system.

### OPTIONAL COMPONENTS, ONLY CHOOSABLE WHEN ORDERING THE FRAME RUNNER

<p><input type="checkbox"/> Saddle abatement system, to facilitate the ascent and descent from the aid.</p> 	<p>Please choose the multiadjustable, padded trunk support size.</p> <p><input type="checkbox"/> M   <input type="checkbox"/> L</p>  <p>In absence of choice, the corresponding size will be supplied.</p>	<p>Please choose the handlebar size.</p> <p><input type="checkbox"/> M   <input type="checkbox"/> L</p>  <p>Please choose the position of the brake on the handlebar or on the ergonomic handles (880).</p> <p><input type="checkbox"/> Left   <input type="checkbox"/> Right</p> <p>In absence of choice, the brake will be mounted on the right.</p>	<p>Please choose the seatpost length. ** The 15 cm (5,9") seatpost does not cover the entire adjustment range.</p> <p><input type="checkbox"/> **15 cm (5,9")</p> <p><input type="checkbox"/> 28 cm (11")</p>  <p>In absence of choice, the 28 cm (11") seatpost will be supplied.</p>
--	---	--	---

Steering limiter spring

Soft

Hard

In absence of choice, size M will be supplied with the soft spring and size L with the hard spring.

### OPTIONAL COMPONENTS, ALSO CHOOSABLE SEPARATELY

#### 101 Wheels dimensions

##### Frame Runner M

In absence of choice, the M size will be supplied with the trio of 24" wheels.

- Front wheel with fork 20" - Rear wheels 24"
- Front wheel with fork 24" - Rear wheels 26"


##### Frame Runner L

In absence of choice, the L size will be supplied with the trio of 26" wheels.

- Front wheel with fork 24" - Rear wheels 24"
- Front wheel with fork 24" - Rear wheels 26"
- Front wheel with fork 26" - Rear wheels 24"
- Front wheel with fork 26" - Rear wheels 28"
- Front wheel with fork 28" - Rear wheels 28"


880 Ergonomic handles adjustable in inclination independently (1 pair). Same size as the Handlebar.

Choose the position of the brake on Ergonomic Handles.




875 Wrist retainer with support and multi-adjustable glove independently.

Left    Right



845 Forearm support.

Left    Right



102 Trunk support wide strap, width adjustable; to be applied to the trunk support.

M    L



876 Monocycle saddle.



103 Monocycle saddle padding.







104 Frame side paddings.



105 Frame rear paddings.



<input type="checkbox"/> <b>106</b> Under saddle padding.	<input type="checkbox"/> <b>107</b> Parking brake, right and left.	<input type="checkbox"/> <b>108</b> Front led lamp. (with USB charging)	<input type="checkbox"/> <b>109</b> Back led lamp. (with USB charging)
<input type="checkbox"/> <b>110</b> Bottle holder. Water bottle not included.	<input type="checkbox"/> <b>111</b> Transparent rear wheels spoke protectors, right and left.	Notes:----- ----- ----- ----- ----- ----- -----	

COMPONENTS MEASURES	<b>Padded trunk support</b>			<b>Handlebar</b>																	
																					
		A	B	C		A															
		<table border="1"> <tr> <td></td> <td style="text-align: center;">A</td> <td style="text-align: center;">B</td> <td style="text-align: center;">C</td> </tr> <tr> <td>M</td> <td style="text-align: center;">16 cm/6,3 in</td> <td style="text-align: center;">32 cm/12,6 in</td> <td style="text-align: center;">5,8 cm/2,3 in</td> </tr> <tr> <td>L</td> <td style="text-align: center;">20 cm/7,9 in</td> <td style="text-align: center;">32 cm/12,6 in</td> <td style="text-align: center;">5,8 cm/2,3 in</td> </tr> </table>		A	B	C	M	16 cm/6,3 in	32 cm/12,6 in	5,8 cm/2,3 in	L	20 cm/7,9 in	32 cm/12,6 in	5,8 cm/2,3 in		<table border="1"> <tr> <td></td> <td style="text-align: center;">A</td> </tr> <tr> <td>M</td> <td style="text-align: center;">55 cm/21,6 in</td> </tr> <tr> <td>L</td> <td style="text-align: center;">64 cm/25,2 in</td> </tr> </table>		A	M	55 cm/21,6 in	L
	A	B	C																		
M	16 cm/6,3 in	32 cm/12,6 in	5,8 cm/2,3 in																		
L	20 cm/7,9 in	32 cm/12,6 in	5,8 cm/2,3 in																		
	A																				
M	55 cm/21,6 in																				
L	64 cm/25,2 in																				
	<b>102</b>			<b>876</b>																	
																					
		A	B		A	B	C														
	<table border="1"> <tr> <td></td> <td style="text-align: center;">A</td> <td style="text-align: center;">B</td> </tr> <tr> <td>M</td> <td style="text-align: center;">8 cm/3,14 in</td> <td style="text-align: center;">45-59 cm/17,7-23,2 in</td> </tr> <tr> <td>L</td> <td style="text-align: center;">8 cm/3,14 in</td> <td style="text-align: center;">63,5-84 cm/25-33 in</td> </tr> </table>		A	B	M	8 cm/3,14 in	45-59 cm/17,7-23,2 in	L	8 cm/3,14 in	63,5-84 cm/25-33 in		<table border="1"> <tr> <td></td> <td style="text-align: center;">A</td> <td style="text-align: center;">B</td> <td style="text-align: center;">C</td> </tr> <tr> <td></td> <td style="text-align: center;">13,5 cm/5,3 in</td> <td style="text-align: center;">31 cm/12,2 in</td> <td style="text-align: center;">9 cm/3,5 in</td> </tr> </table>		A	B	C		13,5 cm/5,3 in	31 cm/12,2 in	9 cm/3,5 in	
	A	B																			
M	8 cm/3,14 in	45-59 cm/17,7-23,2 in																			
L	8 cm/3,14 in	63,5-84 cm/25-33 in																			
	A	B	C																		
	13,5 cm/5,3 in	31 cm/12,2 in	9 cm/3,5 in																		