

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:



Product marked **CE**



**Attention!**  
Read carefully the user and maintenance manual

### FEATURES

- **Ergonomic and modular tricycle** that helps the user to keep a proper posture and makes the ride easier thanks to its design and components.
- **Lowered frame** for easy ascent and descent.
- **Front positioned bottom bracket** for easier pedaling.
- **Pinion** to facilitate the pedaling in safety.
- **Folding frame** to save space during transport.
- **Adjustable handlebar** in height, backward-forward and reversible.
- **Wide and ergonomic handles**, to increase comfort.
- **Bell.**
- **Front and back led lamp.**
- **Front brake with locking device** for ascent and descent, that can be operated by the lever either on the right or on the left side of the handlebar. The distance between the brake lever and the handlebar is adjustable by screw.
- **Adjustable ergonomic saddle** in height, tilt and backward-forward
- **Rear storage basket.**
- **20" pneus wheels.** A mechanism on the rear wheels allows easy adjustment of the tension of the chain.
- **Quick release front wheel.**


Please check off the interested products:  \_\_\_\_\_

### MEASURES

FRAME COLOUR:


- Blue
- Red
- Green

View from above




Storage basket  
Internal measures

L	55,3 cm	in 21,7
M	27,4 cm	in 10,7
N	9,6 cm	in 3,7
Max. load	35 kg	lb 77



C1= C2-2cm



C2= leg length

	Max. Load	Aid Weight	A	B	C1*	D	E	F	G	I	User Height (h)
<b>XL</b>	kg. 100 lb. 220	kg. 28 lb. 61,7	cm 64 in 25,2	cm 82 in 32,2	cm 73-87 in 28,7-34	cm 16 in 6,3	cm 168 in 66	cm 61,5 in 24,2	cm 96 in 37,7	cm 115 in 45,2	cm 155-190 in 61-74,8

\*With standard saddle. You can further adjust the height of the saddle (3cm ±) acting on the rail clamp.

INGLESE - cod.220511 - 06-11-24

Choose the configuration

Fixed Pinion



Free Wheel

Differential and Fixed Pinion













Differential and Free Wheel

Three-speed transmission with differential. Allows braking by counter-pedaling, shifting when the tricycle is stationary, and reversing.

### OPTIONAL COMPONENTS, ONLY CHOOSABLE WHEN ORDERING THE TRICYCLE

<p><input type="checkbox"/> <b>864</b> Disc brake. only for models with differential and without gearbox.</p> 	<p><input type="checkbox"/> <b>847</b> Cranks adjustable in length decreases the distance between saddle and pedal of 3,5 and 5 cm.</p> 	<p><b>Standard handlebaar</b></p> <p><input type="checkbox"/> size L    <input type="checkbox"/> size XL</p> <p>Please choose the position of the brake</p> <p><input type="checkbox"/> Left    <input type="checkbox"/> Right</p> 	<p>In the absence of indications, the brake will be mounted on the right.</p>
---	---	---	---

### OPTIONAL COMPONENTS, ALSO CHOOSABLE SEPARATELY

<p><input type="checkbox"/> <b>880</b> Ergonomic handles adjustable in inclination independently (1 pair). Same size as the Handlebaar.</p> <p>Choose the position of the brake on Ergonomic Handles</p> 	<p><b>845</b> Forearm support. Size L-XL</p> <p><input type="checkbox"/> Left    <input type="checkbox"/> Right</p> 	<p><b>875</b> Wrist retainer with support and multi-adjustable glove independently. Size L-XL</p> <p><input type="checkbox"/> Left    <input type="checkbox"/> Right</p> 	<p><input type="checkbox"/> <b>842P</b> Pelvic support with belt, adjustable in height. Also select the component 883</p> <p><input type="checkbox"/> M    <input type="checkbox"/> L-XL</p> 
<p><b>842T</b> Trunk support, adjustable in height, inclination and backward-forward with belt and joint. Also select the component 883</p> <p><input type="checkbox"/> M    <input type="checkbox"/> L-XL</p> 	<p><input type="checkbox"/> <b>883</b> Rod adjustable in height for Thoracic/Pelvic support (single rod). Please select even the type of support.</p> 	<p><input type="checkbox"/> <b>877</b> Flat backrest</p> 	<p><input type="checkbox"/> <b>843</b> Abductive saddle with holder for supports</p> 
<p><input type="checkbox"/> <b>876</b> Monocycle saddle with holder for supports</p> 	<p><input type="checkbox"/> <b>840</b> Pedal cages with straps (1 pair).</p> 	<p><b>841</b> Pedal cage with length and width adjustable foot imprints and fastening straps, with support kit.</p> <p><input type="checkbox"/> Left    <input type="checkbox"/> Right</p> 	<p><b>846</b> Ankle support adjustable in height, with support kit. Only with component 841.</p> <p><input type="checkbox"/> Left    <input type="checkbox"/> Right</p> <p><input type="checkbox"/> S-M    <input type="checkbox"/> L-XL</p> 
<p><input type="checkbox"/> <b>881</b> Steering limiter.</p> 	<p><input type="checkbox"/> <b>882</b> Steering limiter with springs to return the steering to the straight position.</p> 	<p><input type="checkbox"/> <b>879</b> Front basket.</p> 	<p><input type="checkbox"/> <b>878</b> 878 MTB tires (3 pcs.) in alternative to standard tyres</p> 

COMPONENTS MEASURES

### 842P

	A
mini/s	8,5-14 cm/in 3,3-5,5
M	8,5-21 cm/in 3,3-8,3
L/XL	8,5-23 cm/in 3,3-9

### 842T

	A
mini/s	19-30,5 cm/in 7,4-12
M	19-40 cm/in 7,4-15,7
L/XL	19-46 cm/in 7,4-18,1

### 842P - 842T

	B	C
mini/s	10,5 cm/in 4,1	23 cm/in 9
M	12,5 cm/in 4,9	31 cm/in 12,2
L/XL	17,5 cm/in 6,8	35 cm/in 13,7

### 843

	A	B
mini/s	31 cm/in 12,2	12,5 cm/in 4,9
M/L/XL	40 cm/in 15,7	14 cm/in 5,5

### 876

A	31 cm 12,2"
B	13,5 cm 5,3"
C	9 cm 3,5"

### 841

841 foot imprint size			
	A	B	C
mini/s	from 5 to 7 cm from 1,9 to 2,7 in	from 12,5 to 16,5 cm from 4,9 to 6,5 in	from 6 to 7 cm from 2,3 to 2,7 in
M/L/XL	from 8,2 to 9,4 cm from 3,2 to 3,7 in	from 20 to 22,5 cm from 7,8 to 8,8 in	from 10,4 to 12 cm from 4 to 4,7 in

### Standard handlebar

	A
mini	45 cm / in 17,7
S	50 cm / in 19,6
M	55 cm / in 21,6
L	60 cm / in 23,6
XL	64 cm / in 25,2

Notes:-----  
-----  
-----  
-----  
-----  
-----  
-----  
-----  
-----  
-----  
-----  
-----