

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:




according to the  
ISO 7176-19

Certificate CERAH France  
Marked 

### FEATURES

- Steel folding **frame**. Painted with non-toxic epoxy vanish.
- Removable **upholstery** covered with fire retardant fabric.
- **Rigid backrest** adjustable in tilt on 4 positions, thanks to a simple quick release mechanism.
- Rigid **seat**.
- **Leg rest** adjustable in tilt with quick lock mechanism.
- Tip-up **footrest** with heel pad.
- **Drum brakes**.
- **Wheels**: front wheels 17,5 cm swivel and unpierceable and rear wheels 30 cm unpierceable.

 **Attention!**  
Read carefully the user  
and maintenance manual

Please check off the interested products:




### MEASURES



OBI-4

Frame Colour: **blue**

Upholstery Colour: **green-pink-blue (combined)**



SIZE	Max. Load	Total Weight	A	B	C	D	E	F	G	L	M	H
4	70 kg 154,3 lb	20,1 kg 44,3 lb	66,5 cm 26,1 in	41 cm 16,1 in	124 cm 48,8 in	105 cm 41,3 in	40 cm 15,7 in	76 cm 29,9 in	32-47 cm 12,5-18,5 in	103 cm 40,5 in	51 cm 20 in	150-175 cm 59-68,8 in

ADDITIONAL COMPONENTS

**835** Headrest with  
 parietal suppor



**816** Side supports



**834** Abduction  
 block



**853** Adjustable  
 harness



**894** 45° pelvic belt



**827** Foot straps



**839** Front handle



**837** Tray without  
 hollow



**819** Canopy



**825** Rain cover  
*(to fix at the canopy 819)*



**818** Thermal cover



**858** Net basket



**812** Directional  
 locks for front wheels



**891** Set of 4 tie  
 down hooks  
 for the forward-facing  
 transport of the pushchair in  
 a motor vehicle (private  
 cars, buses..).



Notes:

-----

-----

-----

-----

-----

-----